

# BUECHEL-PAPPAS™

## HIGH FLEX TOTAL KNEE REPLACEMENT SYSTEM

### Primary Knee System



#### Ultracoat - TiN Coated

- Titanium Femur & Tibia
- Lighter, Stronger & Longer lasting

#### High Flex Design



#### Cobalt Chrome

- Cobalt chrome Femur & Tibia

#### Deep Dish Design for Low Contact Strees



#### All Poly

- Medial Pivot model

#### Tibal Pin to prevent spin-out



### Revision Knee System

#### Revision with Wedge

Available in

**Sizes:**  
1, 2, 3, 4 & 5

**Thickness:**  
5 mm, 10 mm &  
15 mm



#### Hinge

Available in

**Sizes:**  
2, 3 & 4



#### Stem Extension



Available in

**Diameter:**  
10 mm, 12 mm,  
14 mm & 16 mm

**Length:**  
50 mm, 75 mm &  
125 mm



**BUECHEL-PAPPAS™ TitaNium knee**  
with  
**Ultracoat®**



**Ultracoat® - the TiN ceramic coating  
with long term history**