



SYSTEM DESCRIPTION - IMPLANT AVAILABILITY & SIZING



Insert Thickness 8, 10, 12, 14, 16 mm



CITIUS[®] Compatibility Chart

| | | | | | | | |
|--|---|-----|-----|---|-----|---|---|
| Femur (L & R) | All Sizes of Femoral Components compatible with all sizes of Tibial Inserts | | | | | | |
| Tibia Sizes | 1 | 2 | 3 | 4 | 5 | 6 | 7 |
| Insert Sizes (Thickness 8, 10, 12, 14, 16 mm) | 1 | | | | | | |
| | | 2/3 | | | | | |
| | | | 4/5 | | | | |
| | | | | | 6/7 | | |
| Patella | All Sizes of Patella Components compatible with all sizes of Femoral Components | | | | | | |

CITIUS



CITIUS[®]

HIGH FLEX TOTAL KNEE REPLACEMENT SYSTEM
POSTERIOR STABILISED







For additional informations kindly contact / scan QR code
TTK Healthcare Ltd.,
 6, Cathedral Road,
 Chennai - 600 086
 P: +91 44 2811 6106 / 08
 E: ortho@ttkhealthcare.com
 W: www.ttkortho.com



www.ttkortho.com

PRODUCT MANUAL



Size distribution according to patients

